


Eukaryotic Topoisomerase I Assay Buffer

Catalog Number:

4010 1ml



Description

•**Topo I Assay Buffer is supplied as a 10x Stock solution.** The buffer is optimized and tested for relaxation and cleavages by topo I. [TopoGEN also offers purified topo enzymes on-line at www.topogen.com.]

Quality Control Analysis

•Vial contains 1ml of 100 mM Tris-HCl, (pH 7.9), 10 mM EDTA, 1.5 M NaCl, 1% BSA, 1 mM spermidine, 50% glycerol. This buffer must be diluted 10 fold for use.

Storage and Shipping Conditions

•The product may be shipped at ambient temperature or on dry ice. Store at 4°C and avoid freezing and thawing as this may inactivate the buffer. At 4°C, this buffer has a useable shelf life of 6 months.

Product Application and Disclaimer

This product is not licensed or approved for administration to humans or animals. It may be used with experimental animals only. The product is for in vitro research diagnostic studies only. The product in non-infectious and non-hazardous to human health. This information is based on present knowledge and does not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship. TopoGEN, Inc. shall not be held liable for product failure due to mishandling and incorrect storage by end user. TopoGEN's liability is limited to credit or product replacement.



TopoGEN, Inc.

For orders: www.topogen.com

For Technical Questions: 614-451-5810


or support@topogen.com

©TopoGEN, Inc. 2009

Eukaryotic Topoisomerase I Assay Buffer

Catalog Number:

4010 1ml



Description

•**Topo I Assay Buffer is supplied as a 10x Stock solution.** The buffer is optimized and tested for relaxation and cleavages by topo I. [TopoGEN also offers purified topo enzymes on-line at www.topogen.com.]

Quality Control Analysis

•Vial contains 1ml of 100 mM Tris-HCl, (pH 7.9), 10 mM EDTA, 1.5 M NaCl, 1% BSA, 1 mM spermidine, 50% glycerol. This buffer must be diluted 10 fold for use.

Storage and Shipping Conditions

•The product may be shipped at ambient temperature or on dry ice. Store at 4°C and avoid freezing and thawing as this may inactivate the buffer. At 4°C, this buffer has a useable shelf life of 6 months.

Product Application and Disclaimer

This product is not licensed or approved for administration to humans or animals. It may be used with experimental animals only. The product is for in vitro research diagnostic studies only. The product in non-infectious and non-hazardous to human health. This information is based on present knowledge and does not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship. TopoGEN, Inc. shall not be held liable for product failure due to mishandling and incorrect storage by end user. TopoGEN's liability is limited to credit or product replacement.



TopoGEN, Inc.

For orders: www.topogen.com

For Technical Questions: 614-451-5810

or support@topogen.com

©TopoGEN, Inc. 2009