

10x Agarose Gel Loading Buffer

Catalog Number:

4038 1ml



Description

- Vial contains 1ml of 10x agarose gel loading buffer (0.005% bromophenol blue, 50% glycerol).

Storage and Shipping Conditions

- The product may be shipped at ambient temperature or on Dry Ice. Store at 4°C.

Product Application and Disclaimer

This product is not licensed or approved for administration to humans or animals. It may be used with experimental animals only. The product is for in vitro research diagnostic studies only. The product is non-infectious and non-hazardous to human health. This information is based on present knowledge and does not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship. TopoGEN, Inc. shall not be held liable for product failure due to mishandling and incorrect storage by end user. TopoGEN's liability is limited to credit or product replacement.



TopoGEN, Inc.

For orders: www.topogen.com

For Technical Questions: 614-451-5810

or support@topogen.com

©TopoGEN, Inc. 2009

10x Agarose Gel Loading Buffer

Catalog Number:

4038 1ml



Description

- Vial contains 1ml of 10x agarose gel loading buffer (0.005% bromophenol blue, 50% glycerol).

Storage and Shipping Conditions

- The product may be shipped at ambient temperature or on Dry Ice. Store at 4°C.

Product Application and Disclaimer

This product is not licensed or approved for administration to humans or animals. It may be used with experimental animals only. The product is for in vitro research diagnostic studies only. The product is non-infectious and non-hazardous to human health. This information is based on present knowledge and does not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship. TopoGEN, Inc. shall not be held liable for product failure due to mishandling and incorrect storage by end user. TopoGEN's liability is limited to credit or product replacement.



TopoGEN, Inc.

For orders: www.topogen.com

For Technical Questions: 614-451-5810

or support@topogen.com

©TopoGEN, Inc. 2009